

DO's for Oxygen Concentrator (OC)



Use OC in well Ventilated Room



Keep OC away from wall/screens



External filters should be cleaned once a week



Use disinfectant/ soap solution to clean OC cabinet



LPM knob should be set at correct amount as prescribed by doctor



Ensure no gap around face mask



Use clean water to fill the humidifier bottle



Follow infection control protocol for disposal of consumables

DON'Ts for Oxygen Concentrator (OC)



Do not use extension cord for power supply



Do not use TAP Water for humidifier bottle



Do not keep OC near curtains

ऑक्सीजन कॉन्सेंट्रेटर (ओ.सी.) के उपयोग के दौरान क्या करें



ओ.सी. का उपयोग हवादार कमरे में करें



ओ.सी. को दीवार/स्क्रीन से दूर रखें



बाहरी फिल्टर को सप्ताह में एक बार साफ करना चाहिए




ओ.सी. कैबिनेट को कीटाणुनाशक / साबुन के घोल से साफ करें



सुनिश्चित करें कि फेस मास्क के आसपास कोई गैप न हो



डॉक्टर द्वारा बताई गई ऑक्सीजन की सही मात्रा पर ही एलपीएम नाॅब को सेट करें



ह्यूमिडिफायर की बोटल भरने के लिए साफ पानी का इस्तेमाल करें



उपयोग में लायी गई सामग्रियों के निपटान के लिए संक्रमण नियंत्रण प्रोटोकॉल का पालन करें

ऑक्सीजन कॉन्सेंट्रेटर(ओ.सी.) के उपयोग के दौरान क्या न करें



बिजली की आपूर्ति के लिए एक्सटेंशन बोर्ड का उपयोग न करें



ह्यूमिडिफायर बोटल के लिए नल के पानी का उपयोग न करें



ओ.सी. को पर्दों के पास न रखें